



一次性使用活体取样钳
**DISPOSABLE BIOPSY
FORCEPS**

www.pengtianmedical.com www.cnpty.com

参数规格 Specification


规格型号 Item No.	钳子杯直径 O.D. (mm)	有效长度 Working Length(mm)	钳口形式 Jaw Shape	是否带针 Spike(N)	是否包塑 Coating (B)	最小可通过通道 Endoscope Channel(mm)
PTY1-23 (16/18/20/23) (T/U/S)(N, B)	2.3	1600/1800/2000/2300	平口Oval(T) 鳄口Alligator(U) 圆杯口Round(S)	是Yes/否No	是Yes/否No	2.8
PTY2-18 (07/10/12/16/18/23) (T/U/S)(N, B)	1.8	700/1000/1200 1600/1800/2300	平口Oval(T) 鳄口Alligator(U) 圆杯口Round(S)	是Yes/否No	是Yes/否No	2.0

Disposable Oval Biopsy Forceps

Item No.	O.D. (mm)	Working Length(mm)	Characteristics	Spike	Coating	Endoscope Channel(mm)
PTBF-23(16/18/20/23)BT	2.3	1600/1800/2000/2300	Oval	No	No	2.8
PTBF-23(16/18/20/23)BN	2.3	1600/1800/2000/2300	Oval	Yes	No	2.8
PTBF-23(16/18/20/23)BTC	2.3	1600/1800/2000/2300	Oval	No	Yes	2.8
PTBF-23(16/18/20/23)BNC	2.3	1600/1800/2000/2300	Oval	Yes	Yes	2.8
PTBF-23(16/18/20/23)BTO	2.3	1600/1800/2000/2300	Oval, distal fine	No	No	2.8
PTBF-23(16/18/20/23)BNO	2.3	1600/1800/2000/2300	Oval, distal fine	Yes	No	2.8
PTBF-23(16/18/20/23)BOC	2.3	1600/1800/2000/2300	Oval, distal fine	No	Yes	2.8
PTBF-23(16/18/20/23)BNCO	2.3	1600/1800/2000/2300	Oval, distal fine	Yes	Yes	2.8

Disposable Alligator Biopsy Forceps

Item No.	O.D. (mm)	Working Length(mm)	Characteristics	Spike	Coating	Endoscope Channel(mm)
PTBF-23(16/18/20/23)BU	2.3	1600/1800/2000/2300	Alligator teeth	No	No	2.8
PTBF-23(16/18/20/23)BUN	2.3	1600/1800/2000/2300	Alligator teeth	Yes	No	2.8
PTBF-23(16/18/20/23)BUC	2.3	1600/1800/2000/2300	Alligator teeth	No	Yes	2.8
PTBF-23(16/18/20/23)BZC	2.3	1600/1800/2000/2300	Alligator teeth	Yes	Yes	2.8
PTBF-23(16/18/20/23)BUO	2.3	1600/1800/2000/2300	Alligator teeth, distal fine	No	No	2.8
PTBF-23(16/18/20/23)BUNO	2.3	1600/1800/2000/2300	Alligator teeth, distal fine	Yes	No	2.8
PTBF-23(16/18/20/23)BUCO	2.3	1600/1800/2000/2300	Alligator teeth, distal fine	No	Yes	2.8
PTBF-23(16/18/20/23)BZCO	2.3	1600/1800/2000/2300	Alligator teeth, distal fine	Yes	Yes	2.8



1.8 series Disposable Biopsy Forceps

Item No.	O.D.(mm)	Working Length(mm)	Characteristics	Spike	Coating	Endoscope Channel(mm)
PTF-18 (07/10/12/16/18/23)BT	1.8	700/1000/1200/1600/1800/2300	Oval	No	No	2.0
PTF-18 (07/10/12/16/18/23)BN	1.8	700/1000/1200/1600/1800/2300	Oval	Yes	No	2.0
PTF-18 (07/10/12/16/18/23)BTC	1.8	700/1000/1200/1600/1800/2300	Oval	No	Yes	2.0
PTF-18 (07/10/12/16/18/23)BNC	1.8	700/1000/1200/1600/1800/2300	Oval	Yes	Yes	2.0
PTF-18 (07/10/12/16/18/23)BU	1.8	700/1000/1200/1600/1800/2300	Alligator teeth	No	No	2.0
PTF-18 (07/10/12/16/18/23)BUN	1.8	700/1000/1200/1600/1800/2300	Alligator teeth	Yes	No	2.0
PTF-18 (07/10/12/16/18/23)BUC	1.8	700/1000/1200/1600/1800/2300	Alligator teeth	No	Yes	2.0
PTF-18 (07/10/12/16/18/23)BZC	1.8	700/1000/1200/1600/1800/2300	Alligator teeth	Yes	Yes	2.0

Disposable Marked Biopsy Forceps

Item No.	O.D.(mm)	Working Length(mm)	Characteristics	Spike	Coating	Endoscope Channel(mm)
PTMF-23(16/18/20/23)BTC	2.3	1600/1800/2000/2300	Oval	No	Yes	2.8
PTMF-23(16/18/20/23)BNC	2.3	1600/1800/2000/2300	Oval	Yes	Yes	2.8
PTMF-23(16/18/20/23)BOC	2.3	1600/1800/2000/2300	Oval, distal fine	No	Yes	2.8
PTMF-23(16/18/20/23)BNCO	2.3	1600/1800/2000/2300	Oval, distal fine	Yes	Yes	2.8
PTMF-23(16/18/20/23)BUC	2.3	1600/1800/2000/2300	Alligator teeth	No	Yes	2.8
PTMF-23(16/18/20/23)BZC	2.3	1600/1800/2000/2300	Alligator teeth	Yes	Yes	2.8
PTMF-23(16/18/20/23)BUCO	2.3	1600/1800/2000/2300	Alligator teeth, distal fine	No	Yes	2.8
PTMF-23(16/18/20/23)BZCO	2.3	1600/1800/2000/2300	Alligator teeth, distal fine	Yes	Yes	2.8
PTMF-18 (07/10/12/16/18/23)BTC	1.8	700/1000/1200/1600/1800/2300	Oval	No	Yes	2.0
PTMF-18 (07/10/12/16/18/23)BNC	1.8	700/1000/1200/1600/1800/2300	Oval	Yes	Yes	2.0
PTMF-18 (07/10/12/16/18/23)BUC	1.8	700/1000/1200/1600/1800/2300	Alligator teeth	No	Yes	2.0
PTMF-18 (07/10/12/16/18/23)BZC	1.8	700/1000/1200/1600/1800/2300	Alligator teeth	Yes	Yes	2.0

Disposable Moving Cups Biopsy Forceps

Item No.	O.D.(mm)	Working Length(mm)	Characteristics	Spike	Coating	Endoscope Channel(mm)
PTCF-23(16/18/20/23)BUC-D	2.3	1600/1800/2000/2300	Alligator teeth, unilateral turn	No	Yes	2.8
PTCF-23(16/18/20/23)BUC-S	2.3	1600/1800/2000/2300	Alligator teeth, bilateral turn	No	Yes	2.8

Reusable Biopsy Forceps 重复性活检钳



Item No. 规格型号	O.D. 钳口直径(mm)	Working Length 有效长度(mm)	Characteristics 特征
PTRF-2316S	2.3	1600	Round cups圆杯口
PTRF-2323S	2.3	2300	Round cups圆杯口
PTRF-2316T	2.3	1600	Oval cups大平口
PTRF-2323T	2.3	2300	Oval cups大平口
PTRF-1816T	1.8	1600	Oval cups平口

产品特点 Product features

- 用于内窥镜下消化道、呼吸道活体组织取样
- This device is used for tissue sampling in the digestive tract or respiratory tract through endoscope.
- The alligator teeth applies for hard tissue sampling and takes an easily tumor tissue.
- With needle it can help accurate sampling.
- Coating outside effectively decreases the friction to the endoscope channel.
- 20cm fine and soft spring at the front end is more suitable for 360 degree turning of endoscope.
- It's ultra-low friction to the endoscope channel can effectively protect the endoscope.
- Moving cups are easy for accurate targeting, no sliding, depressing lesions.



诸暨市鹏天医疗器械有限公司
ZHUJI PENG TIAN MEDICAL INSTRUMENT CO., LTD.

电话Tel: 86-575-87128772 手机Mob: 86-18657572931(INT Market国外销售) 86-13777302827(China国内销售) 传真Fax: 86-575-87128771
网址Web: www.pengtianmedical.com www.cnpty.com 邮箱Email: allen@pengtianmed.com 地址：中国浙江诸暨暨阳经济开发区津津路8号
Add: No.8 Jinjin Road, Jiyang Economic Development District 311800, Zhuji, Zhejiang, China



一次性使用电圈套器

DISPOSABLE POLYPECTOMY
SNARE

www.pengtianmedical.com www.cnpty.com

DISPOSABLE POLYPECTOMY SNARE 一次性使用电圈套器


参数规格 Specification

Item No. 型号	O.D. 外管直径(mm)	Working Length 工作长度(mm)	Loop Diameter 圈径(mm)	Loop Type 圈套器形状	Endoscope Channel 最小可通过钳道(mm)
PTS-18 (12/18) -(15/20/25)-(A/B/C)	1.8	1200/1800	15/20/25	Oval椭圆形(A) Hexagonal六角形(B) Crescent半月形(C)	2.0
PTS-23 (16/18/23)-(10/15/20/25/30/35) -(A/B/C)	2.3	1600/1800/2300	10/15/20/25/30/35	Oval椭圆形(A) Hexagonal六角形(B) Crescent半月形(C)	2.8

产品特点 Product features

Disposable Polypectomy Snare 一次性使用电圈套器

- 与高频电流配合使用，用于经内窥镜对消化道、呼吸道的息肉的切除
This device is used with high frequency to resect small polyp under endoscope.
- 为适应不同息肉形状的需要，有椭圆形、六角形、半月形不同套圈形状
Oval, hexagonal and crescent type of snare are supplied to meet the needs of different shapes of polyp.
- 连接牢固、可旋转操作，操控性好，可实现单手操作。
Connect firmly, available for rotation operation, and well controllability , can be operated by single hand.



诸暨市鹏天医疗器械有限公司
ZHUJI PENG TIAN MEDICAL INSTRUMENT CO., LTD.

电话Tel: 86-575-87128772 手机Mob: 86-18657572931(INT Market国外销售) 86-13777302827(China国内销售) 传真Fax: 86-575-87128771
网址Web: www.pengtianmedical.com www.cnpty.com 邮箱Email: allen@pengtianmed.com 地址：中国浙江诸暨经济开发区津津路8号
Add: No.8 Jinjin Road, Jiyang Economic Development District 311800, Zhuji, Zhejiang, China



一次性使用注射针
DISPOSABLE
INJECTION NEEDLE

www.pengtianmedical.com www.cnpty.com

DISPOSABLE
INJECTION NEEDLE 一次性使用注射针


参数规格 Specification

Item No. 型号	O.D. 外管直径(mm)	Working Length 工作长度(mm)	Needle Diameter 针径(mm)	Needle Length 针长(mm)	Characteristics 特征
PTI-18 (16/18/23)-(22/23/25) (04/05/06)(A/B/C/D)	Φ1.8	1600/1800/2300	22G(0.7) 23G(0.6) 25G(0.5)	4/5/6	Spring Inside & Metal Head (A)有弹簧和金属头(A) Metal Head (B)金属头(B) W/o metal head (C)普通款(C) Coated spring tube (D)包塑弹簧管(D)
PTI-23 (16/18/23)-(22/23/25) (04/05/06) (A/B/C/D)	Φ2.3	1600/1800/2300	22G(0.7) 23G(0.6) 25G(0.5)	4/5/6	Spring Inside & Metal Head (A)有弹簧和金属头(A) Metal Head (B)金属头(B) W/o metal head (C)普通款(C) Coated spring tube (D)包塑弹簧管(D)

产品特点 Product features

Disposable Injection Needle
一次性使用注射针

- 用于经内窥镜对人体腔道进行液体输送、灌洗喷射等操作
This device is used for entering into human body through endoscope to transport liquid to body cavity or lush body cavity or other operations.
- 金属头、弹簧设计减少对工作钳道损害
The metal head and spring reduce the damage to the channel.
- 端面圆角处理，减少对腔道组织的损伤
The round end is designed to reduce damage to the cavity.



诸暨市鹏天医疗器械有限公司
ZHUJI PENG TIAN MEDICAL INSTRUMENT CO., LTD.

电话Tel: 86-575-87128772 手机Mob: 86-18657572931(INT Market国外销售) 86-13777302827(China国内销售) 传真Fax: 86-575-87128771
网址Web: www.pengtianmedical.com www.cnpty.com 邮箱Email: allen@pengtianmed.com 地址：中国浙江诸暨暨阳经济开发区津津路8号
Add: No.8 Jinjin Road, Jiyang Economic Development District 311800, Zhuji, Zhejiang, China



异物爪钳 / 网兜

DISPOSABLE GRASPING FORCEPS
RETRIEVAL NET

www.pengtianmedical.com www.cnpty.com

DISPOSABLE GRASPING
FORCEPS 异物爪钳

参数规格 Specification

Item No. 规格型号	O.D. 外管直径(mm)	Working Length 有效长度(mm)	Prong 爪脚数量
PTG-2316-Z3	2.3	1600	3
PTG-2316-Z4	2.3	1600	4
PTG-2316-Z5	2.3	1600	5
PTG-2318-Z3	2.3	1800	3
PTG-2318-Z4	2.3	1800	4
PTG-2318-Z5	2.3	1800	5
PTG-2323-Z3	2.3	2300	3
PTG-2323-Z4	2.3	2300	4
PTG-2323-Z5	2.3	2300	5

- 用于经内窥镜抓取切下的息肉、病变组织、团状异物
This device is used for grasping the cut polyps, lesion tissue or mass foreign body.

RETRIEVAL NET 网兜

参数规格 Specification

Item No. 规格型号	Outer Diameter 外管直径(mm)	Working Length 有效长度(mm)	Opening Diameter 开幅直径(mm)
PTW-2323-30	2.3	2300	30

- 带网设计，可准确牢固地抓取异物
This device can accurately grasp the foreign body firmly with a net design.



诸暨市鹏天医疗器械有限公司
ZHUJI PENG TIAN MEDICAL INSTRUMENT CO., LTD.

电话Tel: 86-575-87128772 手机Mob: 86-18657572931(INT Market国外销售) 86-13777302827(China国内销售) 传真Fax: 86-575-87128771
网址Web:www.pengtianmedical.com www.cnpty.com 邮箱Email:allen@pengtianmed.com 地址：中国浙江诸暨暨阳经济开发区津津路8号
Add: No.8 Jinjin Road, Jiyang Economic Development District 311800, Zhuji, Zhejiang, China



Holland Medical



一次性使用取石网篮

DISPOSABLE STONE
EXTRACTION BASKET

www.pengtianmedical.com www.cnpty.com

DISPOSABLE STONE EXTRACTION BASKET 一次性使用取石网篮

参数规格 Specification

Item No. 规格型号	O.D. 外管直径(mm)	Working Length 工作长度(mm)	Basket Shape 网篮类型	Basket Width 网篮开幅(mm)	Endoscope Channel 最小可通过钳道(mm)
PTQ-18 (07/10/18) -(20/25/30)-(A/B/C)	1.8	700/1000/1800	20/25/30	Diamond 钻石形(A) Oval 椭圆形(B) Spiral 螺旋形(C)	2.0
PTQ-23 (16/18/23)-(20/25/30)-(A/B/C)	2.3	1600/1800/2300	20/25/30	Diamond 钻石形(A) Oval 椭圆形(B) Spiral 螺旋形(C)	2.8

产品特点 Product features

- 用于抓取内窥镜下结石、息肉切除术后息肉组织或异物
This device is used to grasp the stones, polyps after removal of polyps or foreign bodies under endoscope.
- 旋转手柄可以实现网篮的360°旋转，有效抓取小结石和粉碎性结石
Rotatable handle can make the basket a 360 °rotation, which effectively captures small stones and crushed stones.
- 不锈钢网篮头拉紧定型能力好，套取后不易脱落
Medical-use stainless steel tightens the basket well, which makes it uneasy to fall off after taking.



诸暨市鹏天医疗器械有限公司
ZHUJI PENG TIAN MEDICAL INSTRUMENT CO., LTD.

电话Tel: 86-575-87128772 手机Mob: 86-18657572931(INT Market国外销售) 86-13777302827(China国内销售) 传真Fax: 86-575-87128771
网址Web: www.pengtianmedical.com www.cnpty.com 邮箱Email: allen@pengtianmed.com 地址：中国浙江诸暨暨阳经济开发区津津路8号
Add: No.8 Jinjin Road, Jiyang Economic Development District 311800, Zhuji, Zhejiang, China



一次性止血夹
DISPOSABLE HEMOCLIP

www.pengtianmedical.com www.cnpty.com

DISPOSABLE HEMOCLIP 一次性止血夹

参数规格 Specification

Item No. 规格型号	Opening Size 开口尺寸L1(mm)	Angle of Tooth 齿部角度A(°)	Opening Type 开口类型	O.D. 外管直径(mm)	Working Length 工作长度(mm)	Category 类别
PTZJ-D-135-26-155	L1≥10	135	Standard标准型	2.6	1550	D:Standard III type标准III型
PTZJ-D-135-26-165	L1≥10	135	Standard标准型	2.6	1650	D:Standard III type标准III型
PTZJ-D-135-26-195	L1≥10	135	Standard标准型	2.6	1950	D:Standard III type标准III型
PTZJ-D-135-26-230	L1≥10	135	Standard标准型	2.6	2300	D:Standard III type标准III型
PTZJ-D-135-26-235	L1≥10	135	Standard标准型	2.6	2350	D:Standard III type标准III型
PTZJ-D-135-26-270	L1≥10	135	Standard标准型	2.6	2700	D:Standard III type标准III型

产品特点 Product features

- 用于在内窥镜下夹合消化道内软组织，消化道止血、标记
This device is used for clamping soft tissue, gastrointestinal bleeding, marking under endoscope.
- 可重复开闭，真正实现重复预夹功能
The clips can be opened and closed repeatedly, which truly realize the function of pre-clamping.
- 可旋转功能,有助于医生以最合适的角度夹持组织
The clips can be rotated to help the doctor hold the tissue at the most appropriate angle.



诸暨市鹏天医疗器械有限公司
ZHUJI PENG TIAN MEDICAL INSTRUMENT CO., LTD.

电话Tel: 86-575-87128772 手机Mob: 86-18657572931(INT Market国外销售) 86-13777302827(China国内销售) 传真Fax: 86-575-87128771
网址Web:www.pengtianmedical.com www.cnpty.com 邮箱Email:allen@pengtianmed.com 地址：中国浙江诸暨暨阳经济开发区津津路8号
Add: No.8 Jinjin Road, Jiyang Economic Development District 311800, Zhuji, Zhejiang, China



一次性清洗刷/重复性清洗刷

DISPOSABLE DOUBLE ENDED
CLEANING BRUSHES
REUSABLE CLEANING BRUSH

www.pengtianmedical.com www.cnpty.com

DISPOSABLE DOUBLE ENDED CLEANING BRUSHES 一次性清洗刷

参数规格 Specification

Item No. 规格型号	O.D. 外管直径(mm)	Working Length 有效长度(mm)	Brush Diameter 刷头直径(mm)
PTCB-A-03-16	1.8	1600	3
PTCB-A-03-18	1.8	1800	3
PTCB-A-03-23	1.8	2300	3
PTCB-A-06-16	1.8	1600	6
PTCB-A-06-18	1.8	1800	6
PTCB-A-06-23	1.8	2300	6
PTCB-A-06/11-16	1.8	1600	6/11
PTCB-A-06/11-18	1.8	1800	6/11
PTCB-A-06/11-23	1.8	2300	6/11

REUSABLE CLEANING BRUSH 重复性清洗刷

参数规格 Specification

Item No. 规格型号	O.D. 外管直径(mm)	Working Length 有效长度(mm)	Brush Diameter 刷头直径(mm)
PTCB-06-12	1.8	1200	6
PTCB-06-16	1.8	1600	6
PTCB-06-18	1.8	1800	6
PTCB-06-23	1.8	2300	6

产品特点 Product features

- 用于清洗内窥镜钳道
Cleaning brushes are used to clean the endoscope channel.
- 一次性双头清洗刷使用方便，高效清洗
Disposable double ended cleaning brushes are used conveniently and effectively.
- 一次性、重复性可供不同选择
Disposable type and reusable type are available for different choices.




诸暨市鹏天医疗器械有限公司
ZHUJI PENG TIAN MEDICAL INSTRUMENT CO., LTD.

电话Tel: 86-575-87128772 手机Mob: 86-18657572931(INT Market国外销售) 86-13777302827(China国内销售) 传真Fax: 86-575-87128771
网址Web: www.pengtianmedical.com www.cnpty.com 邮箱Email: allen@pengtianmed.com 地址：中国浙江诸暨暨阳经济开发区津津路8号
Add: No.8 Jinjin Road, Jiyang Economic Development District 311800, Zhuji, Zhejiang, China



Holland
Medical



一次性使用细胞刷/水汽道刷
DISPOSABLE CYTOLOGY BRUSHES
PORT HOLE CLEANING BRUSH (B)

www.pengtianmedical.com www.cnpty.com

DISPOSABLE CYTOLOGY BRUSHES

一次性使用细胞刷

参数规格 Specification

Item No. 规格型号	O.D. 外管直径(mm)	Working Length 有效长度(mm)	Brush Head Model 刷头型号	Brush Diameter 刷头直径(mm)	Bristle Length 刷长(mm)	Endoscope Channel 最小可通过钳道(mm)
PTX-18 (10/12/16/18) -I- (20/30/40/50)10	1.8	1000/1200/1600/1800	Straight直型刷头(I)	2/3/4/5	10	2.0
PTX-23 (16/18/23) -I- (20/30/40/50)10	2.3	1600/1800/2300	Straight直型刷头(I)	2/3/4/5	10	2.8

- 用于内窥镜下消化道呼吸道组织取样
The device is used for endoscopic tissue sampling of the digestive tract or respiratory tract.
- 球形光滑头端，有效降低对腔道壁的刺激和损伤
Spherical smooth head effectively reduces the stimulation and damage of the cavity wall.
- 直形细胞刷使得容易进入消化道呼吸道深处
Straight type makes it easy to enter the depth of the digestive tract or respiratory tract.

PORT HOLE CLEANING BRUSH (B) 水汽道刷

参数规格 Specification

Item No. 规格型号	Brush Diameter 刷头直径(mm)
PTCB-B	5/11

- 用于内镜水气按钮管道清洗
The device is used for cleaning endoscopic water gas button pipe.



诸暨市鹏天医疗器械有限公司
ZHUJI PENG TIAN MEDICAL INSTRUMENT CO., LTD.

电话Tel: 86-575-87128772 手机Mob: 86-18657572931(INT Market国外销售) 86-13777302827(China国内销售) 传真Fax: 86-575-87128771
网址Web: www.pengtianmedical.com www.cnpty.com 邮箱Email: allen@pengtianmed.com 地址: 中国浙江诸暨暨阳经济开发区津津路8号
Add: No.8 Jinjin Road, Jiyang Economic Development District 311800, Zhuji, Zhejiang, China



Holland Medical



活检阀 BIOPSY VALVE

使内窥镜附件顺畅的进入钳道，并保护钳道，减少内窥镜损伤
This device is used to make endoscopic accessory smoothly enter the clamp and avoid injury of endoscope channel.

www.pengtianmedical.com www.cnpty.com

BITE BLOCK 口垫



PTB-YO-01



PTB-YO-02



PTB-YO-03



PTB-YO-04



PTB-YO-05



PTB-YO-06

Bite Block口垫

- 表面光滑，有效保护患者口腔和内窥镜钳道不受损伤
Smooth surface can effectively protect oral cavity of patients and avoid damage to endoscope channel.
- 可选配松紧带，增加安全性和便利度
An optional elastic cord is included for added safety and convenience.
- 不同规格串适合大人小孩的不同需求
Different sizes suit different needs of adults and children.



诸暨市鹏天医疗器械有限公司
ZHUJI PENG TIAN MEDICAL INSTRUMENT CO., LTD.

电话Tel: 86-575-87128772 手机Mob: 86-18657572931(INT Market国外销售) 86-13777302827(China国内销售) 传真Fax: 86-575-87128771
网址Web: www.pengtianmedical.com www.cnpty.com 邮箱Email: allen@pengtianmed.com 地址: 中国浙江诸暨暨阳经济开发区津津路8号
Add: No.8 Jinjin Road, Jiyang Economic Development District 311800, Zhuji, Zhejiang, China