



ISTERIX



Thermoplastic rubber oval balloon of 2,5 cc on the distal tip to ensure the positioning of the catheter

Placement sheath



Shapeable HS CATHETER

Hysterosalpingography balloon catheter

- Placement sheath decreases the need for a tenaculum
- Allows for clear view of the external cervical os
- Stylet available for cases when additional stiffness is needed

CODES TO ORDER

CODE	DIAMETER (F)	CATHETER	PACKAGE
VMI1154	5	standard with balloon	10
VMI1155	7	standard with balloon	10
VMI1156	NO	stylet	10
VMI1185	5	shapeable with balloon	10
VMI1187	7	shapeable with balloon	10
VMI1188	5	front delivery	10
VMI1189	7	front delivery	10
VMI2002	5	shapeable sheath	10